



AlphaID™ At Home

Genetic Health Risk Service

Package Insert



: PROGENIKA BIOPHARMA S.A.
Parque Tecnológico de Bizkaia
Ibaizabal bidea, Edificio 504
48160 Derio – Bizkaia - SPAIN

GRIFOLS

Contents:

- [Intended Use](#)
- [Consideration for testing](#)
- [Summary and explanation of the Service](#)
- [Important](#)
- [Warnings and limitations of the Service](#)
- [For healthcare professionals](#)
- [Clinical Performance](#)
- [Analytical Performance](#)
- [User studies](#)
- [References](#)
- [Trademarks](#)
- [Customer Support](#)

Intended use:

The Alphalid™ At Home Genetic Health Risk Service uses qualitative genotyping to detect clinically relevant genetic variants associated with alpha-1 antitrypsin deficiency (AATD) in genomic DNA isolated from human saliva collected from individuals ≥ 18 years with ORAcollect®·Dx (OCD-100.014) for the purpose of reporting and interpreting Genetic Health Risks (GHR).

This Service is indicated for reporting 14 genetic variants in the *SERPINA1* gene: PI*S; PI*Z; PI*I; PI*M procida; PI*M malton; PI*S iiyama; PI*Q0 granite falls; PI*Q0 west; PI*Q0 bellingham; PI*F; PI*P lowell; PI*Q0 mattawa; PI*Q0 clayton, and PI*M heerlen. The report describes if a person is at an increased risk of developing either lung and/or liver disease linked to AATD. The Service does not describe a person's overall risk of developing lung and/or liver disease. AATD is more common in persons of European descent.

Consideration for testing:

- The COPD Foundation, the World Health Organization, the American Thoracic Society, and the GOLD COPD guidelines recommend testing all COPD patients for variants linked to AATD ([1](#), [2](#), [3](#), [4](#)). Testing also extends to people with unexplained liver disease ([1](#), [3](#)).
- Signs and Symptoms Linked to AATD includes:
 - Shortness of breath and wheezing
 - Chronic cough
 - Lung disease, including chronic obstructive pulmonary disease (COPD), emphysema, chronic bronchitis, and bronchiectasis
 - Liver disease, including cirrhosis, jaundice, hepatic enzyme elevations, chronic hepatitis, and liver scarring (fibrosis)
- To find a doctor with experience in AATD from the Alpha-1 Foundation's Clinical Resource Center, visit <http://www.alpha1.org/alphas-friends-family/resources/find-an-alpha-1-specialist>

- To find a healthcare professional near you with experience testing for AATD, please visit <http://www.AlphaFindADoctor.com>
- To find a genetic counselor, visit <https://findageneticcounselor.nsgc.org>

Summary and explanation of the Service:

The AlphaID™ At Home Genetic Health Risk Service is composed of:

- AlphaID™ At Home Saliva Collection Kit: for human saliva sample collection with ORAcollect®·Dx (OCD-100.014).
- A1AT Genotyping Test: for the genetic analysis and detection of genetic variants associated with AATD.
- AlphaID™ At Home Genetic Health Risk Service website and results portal (www.AlphaIDAtHome.com): to provide the contents and the procedure to order and use the Service.

The consumers order a kit and collect their sample at home. After registering the kit in a HIPAA compliant website portal, they mail their sample to a third-party laboratory in the pre-paid shipping box. When their result is ready, they will receive an email telling them to log onto their secure account to review their results.

Important:

- Please follow the instructions in the AlphaID™ At Home Saliva Collection Kit for sample collection with ORAcollect®·Dx (OCD-100.014) and use the collection device within the expiration date to ensure your DNA results can be processed and connected to your online account.
- If you have a family history of AATD, talk to a healthcare professional about family testing.

Warnings and limitations of the Service:

- DTC – Direct-To-Consumer.
- It is intended for users ≥18 years old.
- It is not a substitute for an appointment with a healthcare professional. We strongly recommend you consult with a healthcare professional if you have any questions or concerns about your result or health.
- It does not diagnose any disease or condition. Only a healthcare professional can diagnose a disease or condition.
- It does not determine if you have or will develop lung and/or liver disease linked to AATD during your lifetime.
- It cannot be used to make healthcare decisions. It does not tell you anything about your overall health. Only a healthcare professional can help you with healthcare decisions.
- It detects 14 variants in the *SERPINA1* gene linked to AATD. These 14 variants explain 95% of AATD cases. It does not detect all possible variants linked to AATD.

- There may be other, non-genetic factors that affect your risk. The Service does not determine your overall risk of developing lung and/or liver disease.
- Different companies offering a genetic risk test may be measuring different genetic variants for the same condition, so you may get different results from a different test.
- The test may not be able to determine a result for all variants analyzed.
- Some people feel a little anxious about getting genetic health results. This is normal. If you feel very anxious, you should speak to your doctor or a genetic counselor prior to collecting your sample for testing. You may also consider getting your test done by your doctor.
- The laboratory may be unable to process every person's sample. The probability that the laboratory cannot process a sample can be up to 0.8%. If this happens, you will receive an email notification. You will also receive another AlphaID™ At Home Saliva Collection Kit to provide a new sample with ORAcollect®·Dx (OCD-100.014) to the laboratory.
- Share your results report with a healthcare professional. Healthcare professionals can answer questions you may have about your results, risk, and how they may apply to your health. They should know you were tested for AATD. You should also inform a healthcare professional if you:
 - Have symptoms of lung or liver disease
 - Have a personal or family history of lung or liver disease
 - Are feeling anxious, uncertain, or concerned about your genetic result or risk
 - Have questions about any risk factors
- The user's race, ethnicity, age and sex may affect how the genetic results are interpreted.

For healthcare professionals:

- This test is not intended to diagnose a disease, determine medical treatment, or tell the user anything about their current state of health.
- This test is intended to provide users with their genetic information to inform health related lifestyle decisions and conversations with healthcare professionals.
- Healthcare professionals should base any diagnostic or treatment decisions on testing and/or other information determined to be appropriate for the patient.

Clinical performance:

Alpha-1 antitrypsin deficiency (AATD) is an under-recognized hereditary disorder. It is passed on from parents to their children through genes and is associated with the premature onset of lung and/or liver disease ([5](#)). More than 90% of the estimated 100,000 people in the United States with AATD don't know they have it ([6](#), [7](#)).

AATD occurs in people of all ethnicities worldwide ([5](#)). However, it is most common in people of European descent. AATD affects about 1 in 1,500 to 3,500 people of European descent ([1](#), [8](#))

Many individuals with AATD are likely undiagnosed, particularly people with a lung condition called chronic obstructive pulmonary disease (COPD). COPD can be caused by AATD; however, AATD is rarely diagnosed. Some people with AATD are misdiagnosed with asthma ([9](#)). Two to

three percent of patients with COPD in the United States are estimated to have AATD ([1](#)). The percentage of patients with liver disease who have AATD is not known.

AATD is caused by certain genetic variants in the *SERPINA1* gene. These variants cause low levels of a protein called alpha-1 antitrypsin (AAT). Low levels of the AAT protein can lead to AATD ([1](#), [2](#)).

A1AT Genotyping Test detects 14 variants in the *SERPINA1* gene linked to AATD: PI*S; PI*Z; PI*I; PI*M procida; PI*M malton; PI*S iiyama; PI*Q0 granite falls; PI*Q0 west; PI*Q0 bellingham; PI*F; PI*P lowell; PI*Q0 mattawa; PI*Q0 clayton and PI*M heerlen. These 14 variants explain 95% of AATD cases ([5](#)). These variants are mainly found in European population, except for PI*S iiyama, which was described in the Asian population ([10](#)).

The most common variants are PI*S and PI*Z. Published studies estimate frequency ranges of 5-10% and 1-3% for PI*S and PI*Z, respectively, in the European population ([8](#)). An analysis of the prevalence of PI*S and PI*Z amongst the five major ethnic subgroups in the United States has demonstrated that the highest risk for AATD is found in Caucasians, followed by Hispanics and Blacks, with the lowest prevalence amongst Mexican Americans and no risk amongst Asian ([8](#)). The remaining 12 genetic variants tested are reported with a very low frequency in the population ([8](#)).

Detailed description of the 14 variants detected by the AlphaID™ At Home Genetic Health Risk Service in the *SERPINA1* gene and the supported References can be found [here](#).

The Service interpretation provides information about the lung and liver disease risk. Detailed description on the risk categories used by the Service can be found [here](#). The risk categorization is done based on the reported clinical cases for each genetic result. References that support the lung and liver disease risk categorization of each genetic result can be found in [Table 1](#).

Analytical performance:

Saliva samples are processed using A1AT Genotyping Test. This test is based on a labeled multiplex PCR amplification and hybridization (Thermal Cycler, Thermo Fisher Scientific) and signal detection using Luminex's xMAP® Technology. The results were analyzed using a proprietary software to obtain the result report file. The result report file is processed by AlphaID™ At Home website and results portal to provide the final report to the user.

Accuracy

A method comparison study was performed to assess the accuracy of A1AT Genotyping Test to correctly detect the genetic variants. A total of 227 samples representing all genetic variants interrogated by the assay were analyzed and compared with Bi-Directional-Sequencing (BDS) (reference method). Percent Agreement (PA) between the two methods for the overall variants and samples (14 variants and 227 samples) was 100% (227/227) with a 95% confidence interval 98.3% to 100%. The percentage of overall "Invalid Tests" was 0% (0/227) with 95% confidence interval 0% to 1.7%.

Precision/Reproducibility

A precision study was performed to assess the reproducibility of the test under different conditions. A total of 792 sample replicates were processed across three sites. The study included two operators per site and was conducted over three days. The reproducibility obtained for the test was 100%.

Minimum DNA input

The minimum DNA input was obtained by testing five human cell line samples using two lots of reagents. The lowest limit of detection was determined to be 0.0215 ng/ μ l.

Interfering substances

Studies were performed to evaluate the potential interference of substances that may be present in saliva samples.

Endogenous substances: Saliva samples collected in ORAcollect®·Dx (OCD-100.014) from five (5) random donors were spiked with substances usually found in saliva samples: salivary α-amylase, hemoglobin, immunoglobin A and total protein. These substances did not interfere with test performance.

Exogenous substances: Saliva samples from five (5) random donors were collected in ORAcollect®·Dx (OCD-100.014) prior to the activity, immediately after the activity and 30 minutes after the activity. Seven (7) activity groups were tested: eating food without beef, eating food with beef, drinking, smoking, chewing gum, using mouthwash and brushing teeth. The results indicated that saliva samples collected in ORAcollect®·Dx (OCD-100.014) should be collected at least 30 minutes after the activity, which is compatible with the instructions for use of ORAcollect®·Dx (OCD-100.014) included in the AlphalD™ At Home Saliva Collection Kit: "Do not eat, drink, smoke or chew gum for 30 minutes before collecting saliva sample."

Microbial substances: Five microbial interfering substances usually found in saliva samples were tested: *Staphylococcus epidermidis*, *Streptococcus mutans*, *Lactobacillus casei*, *Actinomyces viscosus* and *Candida albicans*. These substances did not interfere with test performance.

Interfering variants

The performance of this test may be affected by the presence of rare variants, such as:

Not tested interfering variants:

rs149537225 for PI*S (rs17580); rs551595739 for PI*Z (rs28929474); rs199422213 for P*I (rs28931570), PI*M procida (rs28931569), PI*M malton (rs775982338) and PI*S iiyama (rs55819880); rs544632177 and rs577164283 for PI*F (rs28929470); rs1049800 for PI*P lowell (rs121912714); rs61761869 and rs372571769 for PI*M heerlen (rs199422209); rs148207011 for PI*Q0 granite falls (rs267606950) and PI*Q0 west (rs751235320); rs200634040 and rs72552401 for PI*Q0 bellingham (rs199422211); rs148362959 and rs372571769 for PI*Q0 mattawa (rs763023697), and rs143329723, rs121912712 and rs372571769 for PI*Q0 clayton (rs764325655).

Tested interfering variants:

- When PI*M malton (c.226_228delTTC, rs775982338) is delTTC/delTTC, the reported result for PI*S iiyama (c.230C>T, rs55819880) can be either a homozygous -/-, a heterozygous +/-, or not detected. In any case, the Sample Result PI*M malton/PI*M malton will be reported.
- For a sample compound heterozygous for PI*M malton (c.226_228delTTC, rs775982338) and PI*S iiyama (c.230C>T, rs55819880), a homozygous +/+ result for both PI*S iiyama and PI*M malton will be provided. The Sample Result PI*M malton/PI*S iiyama will be reported. Either sample result (PI*M malton/PI*S iiyama; or PI*M malton/PI*M malton; or PI*S iiyama/PI*S iiyama) is phenotypically associated with a severe A1AT deficiency.
- Variant PI*Q0 amersfoort (c.552C>G, rs19942210) (reference 11) in homozygosity affects to the detection of variant PI*Q0 granite falls (c.552delC, rs267606950) leading to a false positive detection (+/+ or +/-) or an Invalid Test sample result. Variant PI*Q0 amersfoort is described with a global minor frequency (MAF) <0.001%. Both PI*Q0 amersfoort/PI*Q0 amersfoort and PI*Q0 granite falls/PI*Q0 granite falls allelic variant combinations are phenotypically associated with a severe A1AT deficiency.
- The intronic variants rs375637084 (c.1066-26C>T) and rs372571769 (c.1066-25G>A) affect to the detection of variant PI*Z, leading to a false +/+ homozygous result (PI*Z/PI*Z) instead of +/- heterozygous result (PI*M/PI*Z). These variants are described with a global minor frequency (MAF) <0.001% and 0.16%, respectively.
- The variant PI*P Donauwoerth (c.1093G>A, rs143370956) (reference 12) affects to the detection of variant PI*Z, leading to a false -/- homozygous result (PI*M/PI*M) instead of +/- heterozygous result (PI*M/PI*Z). This variant is described with a global minor frequency (MAF) 0.06%.
- The variant rs1049800 (c.840T>C) affects to the detection of variant PI*P lowell, leading to a false +/+ homozygous result (PI*P lowell/PI*P lowell) instead of +/- heterozygous result (PI*M/PI*P lowell). This variant is described with a global minor frequency (MAF) 5.7%.
- Variant PI*Q0 bolton (c.1158delC, rs764325655) (reference 13) in homozygosity affects to the detection of variant PI*Q0 clayton (c.1158dupC, rs764325655) leading to a false positive detection (+/+ or +/-) or an Invalid Test sample result. Variant PI*Q0 bolton is described with a global minor frequency (MAF) <0.001%. Both PI*Q0 bolton/PI*Q0 bolton and PI*Q0 clayton/PI*Q0 clayton allelic variant combinations are phenotypically associated with a severe A1AT deficiency.
- The intronic variant rs1019260714 (c.917+55T>C) affects to the detection of variant PI*S, leading to a false +/+ homozygous result (PI*S/PI*S) instead of +/- heterozygous result (PI*M/PI*S). This variant is described with a global minor frequency (MAF) <0.001%.
- The variant c.1178C>G (rs199422209) in homozygosity affects to the detection of the variant PI*M heerlen (c.1178C>T, rs199422209), leading to a false positive detection (+/-) or an Invalid Test sample result. This variant is described with a global minor frequency (MAF) <0.001%.
- The variant rs201774333 (c.1095C>G) affects to the detection of variant PI*Z, leading to a false +/+ homozygous result (PI*Z/PI*Z) instead of +/- heterozygous result (PI*M/PI*Z). This variant is described with a global minor frequency (MAF) 0.04%.

For more information of the variants, introduce the rs ID at: <https://www.ncbi.nlm.nih.gov>

User studies:

AlphaID™ At Home Saliva Collection Kit user study for ORAcollect®-Dx (OCD-100.014)

Saliva collection kit user study was performed to assess how users understand the saliva collection kit instructions and to assess the ability of naïve users to provide samples adequate for the Service. Study participants represented a demographically diverse US population of naïve users (389 participants). Participants collected and mailed a saliva sample and answered a survey about the collection kit instructions. Saliva samples were processed according to A1AT Genotyping Test package insert. The overall comprehension rate on the collection kit instructions was 98.6%. A1AT Genotyping Test result was obtained for 386 of 389 samples (99.2%) Therefore, the probability that the laboratory cannot process a sample can be up to 0.8%.

AlphaID™ At Home Genetic Health Risk Service report user comprehension study

The user comprehension study for the AlphaID™ At Home Genetic Health Risk Service showed that a demographically diverse US population of naïve users (525 participants) of the Service reports had excellent comprehension of the service's purpose, limitations, results, relevance of ethnicity, other factors that may impact test results, and appropriate next steps. Comprehension was tested through a two-step process. First, participants' comprehension was tested prior to viewing the educational module and Service reports. Second, participants were shown the educational module and the Service reports. Participants completed a survey after the first and second step. As a result, each comprehension domain achieved a minimum of 90.1% or higher user comprehension score in the first step, and 94.0% or higher user comprehension score in the second step, across all reports. The overall comprehension scores were of 92.7% and 96.8% across all comprehension domains and reports, for the first and second step respectively.

References:

General references

1. American Thoracic Society; European Respiratory Society. American Thoracic Society/European Respiratory Society statement: standards for the diagnosis and management of individuals with alpha-1 antitrypsin deficiency. Am J Respir Crit Care Med. 2003 Oct 1;168(7):818-900.
2. World Health Organization. Alpha 1-antitrypsin deficiency: memorandum from a WHO meeting. Bull World Health Organ. 1997;75(5):397-415.
3. Sandhaus RA, Turino G, Brantly ML, Campos M, Cross CE, Goodman K, Hogarth DK, Knight SL, Stocks JM, Stoller JK, Strange C, Teckman J. The Diagnosis and Management of Alpha-1 Antitrypsin Deficiency in the Adult. Chronic Obstr Pulm Dis. 2016 Jun 6;3(3):668-682.
4. Global strategy for the diagnosis, management, and prevention of Chronic Obstructive Pulmonary Disease. GOLD COPD 2021 report.

5. Stoller JK, Hupertz V, Aboussouan LS. Alpha-1 Antitrypsin Deficiency. Deficiency. 2006 Oct 27 [Updated 2020 May 21]. In: Adam MP, Ardinger HH, Pagon RA, et al., editors. GeneReviews [Internet]. Seattle (WA): University of Washington, Seattle; 1993-2020.
6. Campos MA, Wanner A, Zhang G, Sandhaus RA. Trends in the diagnosis of symptomatic patients with alpha1-antitrypsin deficiency between 1968 and 2003. Chest. 2005 Sep;128(3):1179-86.
7. Chorostowska-Wynimko J. Targeted screening programmes in COPD: how to identify individuals with α1-antitrypsin deficiency. Eur Respir Rev. 2015 Mar;24(135):40-5.
8. de Serres F, Blanco I. Role of alpha-1 antitrypsin in human health and disease. J Intern Med. 2014 Oct;276(4):311-35.
9. Genetics Home Reference-NIH (<https://ghr.nlm.nih.gov/condition/alpha-1-antitrypsin-deficiency#statistics>)
10. Seyama K. State of alpha1-antitrypsin deficiency in Japan. Respirology. 2001 Jun;6 Suppl:S35-8.
11. Prins J, van der Meijden BB, Kraaijenhagen RJ, Wielders JP. Inherited chronic obstructive pulmonary disease: new selective-sequencing workup for alpha1-antitrypsin deficiency identifies 2 previously unidentified null alleles. Clin Chem. 2008 Jan;54(1):101-7.
12. Faber J.P, Poller W, Weidinger S, Kirchgesser M, Schwaab R, Bidlingmaier F, Olek K. Identification and DNA sequence analysis of 15 new alpha1- antitrypsin variants, including two PI*QO alleles and one deficient PI*M allele. Am J Hum Genet. 1994. 55:1113-21.
13. Fraizer GC, Siewertsen M, Harrold TR, Cox DW. Deletion/frameshift mutation in the alpha1-antitrypsin null allele, PI*QObolton. Hum Genet. 1989 Nov;83(4):377-82.

References that support the lung and liver disease risk categorization

Table 1. References that support the lung and liver disease risk categorization for each possible genetic result of the AlphalD™ At Home Genetic Health Risk Service. The most frequent and/or most studied genetic results (No Variants, PI*S, PI*Z, PI*S and PI*S, PI*S and PI*Z, PI*Z and PI*Z) are showed in the first rows of the table. The complete list of bibliographic references can be found at the bottom of the table.

Genetic Result in the AlphalD™ At Home Genetic Health Risk Service Report	Number of variants	Reported Lung Disease Risk Category	References that support Lung Disease Risk Category	Reported Liver Disease Risk Category	References that support Liver Disease Risk Category
No Variants	0	Not likely at risk for AATD	ATS/ERS 2003 ; DeMeo & Silverman 2004 ; de Serres & Blanco 2014 ; Stoller et al 2006 (updated 2017) ; Stoller et al 2006 (updated 2020)	Not likely at risk for AATD	ATS/ERS 2003 ; de Serres & Blanco 2014
PI*S	1	Not likely at increased risk ¹	Matzen et al 1977 ; Chang-Yeung et al 1978 ; Ostrow et al 1978 ; Gulsvik et al 1979 ; Dahl et al 2002	Not likely at increased risk	Eigenbrodt et al 1997 ; Graziadei et al 1998 ; Arnaud et al 1977 ; Bell et al 1990 ; Carlson et al 1981

Genetic Result in the AlphaID™ At Home Genetic Health Risk Service Report	Number of variants	Reported Lung Disease Risk Category	References that support Lung Disease Risk Category	Reported Liver Disease Risk Category	References that support Liver Disease Risk Category
PI*Z	1	Not likely at increased risk ¹ (never-smokers)	Seersholt et al 2000 ; Klayton et al 1975 ; Matzen et al 1977 ; Chang-Yeung et al 1978 ; Gulsvik et al 1979 ; Dahl et al 2002	Slightly increased risk ³	Schneider et al 2020 ; Ferrarotti et al 2005 ; Bell et al 1990 ; Eigenbrodt et al 1997 ; Graziadei et al 1998 ; Blenkinsopp & Haffenden 1977 ; Strnad et al 2019 ; Schaefer et al 2018 ; Propst et al 1992 ; Arnaud et al 1977
		Slightly increased risk (ever-smokers)	Sørheim et al 2010 ; Molloy et al 2014		
PI*S and PI*S	2	Not likely at increased risk	Dahl et al 2002 ; Fagerhol et al 1969 ; Gulsvik et al 1979 ; Kueppers et al 1977 ; Lieberman et al 1986 ; Lochon et al 1978 ; Sandford et al 1999 ; Arnaud et al 1977	Not likely at increased risk	Eigenbrodt et al 1997 ; Arnaud et al 1977 ; Strnad et al 2019 ; Carlson et al 1981
PI*S and PI*Z	2	Slightly increased risk	McElvaney et al 2020 ; Fagerhol & Hauge 1969 ; Abboud et al 1979 ; Bartmann et al 1985 ; Lieberman et al 1986 ; Gulsvik et al 1979 ; Dahl et al 2002	Slightly increased risk ³	McElvaney et al 2020 ; Ferrarotti et al 2005 ; Eigenbrodt et al 1997 ; Graziadei et al 1998 ; Propst et al 1992 ; Strnad et al 2019
PI*Z and PI*Z	2	Increased risk	McElvaney et al 2020	Slightly increased risk ³	McElvaney et al 2020 ; Ferrarotti et al 2005 ; Schneider et al 2020 ; Hamesch et al 2019 ; Eriksson et al 1986 ; Elzouki & Eriksson 1996 ; Graziadei et al 1998 ; Propst et al 1992 ; Bell et al 1990
PI*I	1	Not likely at increased risk ¹	Arnaud et al 1978 ; Baur & Bencze 1987 ; Duk et al 2016 ; Zhumagaliev et al 2017	Not likely at increased risk	Arnaud et al 1978 ; Baur & Bencze 1987
PI*M procida	1	Not likely at increased risk ¹	Montealegre et al 2006	Not likely at increased risk ⁴	Ferrarotti et al 2005 ; Balduyck et al 2014

Genetic Result in the AlphaID™ At Home Genetic Health Risk Service Report	Number of variants	Reported Lung Disease Risk Category	References that support Lung Disease Risk Category	Reported Liver Disease Risk Category	References that support Liver Disease Risk Category
PI*M malton	1	Not likely at increased risk ¹	Sproule et al 1983 ; Allen et al 1986 ; Canva et al 2001 ; Orrù et al 2005 ; Corda et al 2006 ; Joly et al 2015 ; Figueira Gonçalves et al 2017	Slightly increased risk ³	Canva et al 2001 ; Janciauskiene et al 2004 ; Orrù et al 2005 ; Corda et al 2006 ; Joly et al 2015 ; Figueira Gonçalves et al 2017 ; Callea et al 2018
PI*S iiyama	1	Not likely at increased risk ¹	Takabe et al 1992	Unknown risk	Clinical cases not reported
PI*Q0 granite falls	1	Not likely at increased risk ¹	Poller et al 1999	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 west	1	Not likely at increased risk ¹	Laubach et al 1993	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 bellingham	1	Not likely at increased risk ¹	Cook et al 1994	Not likely at increased risk ⁴	Clinical cases not reported
PI*F	1	Not likely at increased risk ¹	Cockcroft et al 1981 ; Beckman et al 1984 ; Cook et al 1996 ; Kwok et al 2004 ; Tete-Benissan & Gbeassor 2011 ; Duk et al 2016	Not likely at increased risk ⁴	Kelly et al 1989 ; Tete-Benissan & Gbeassor 2011
PI*P lowell	1	Not likely at increased risk ¹	Seri et al 1992 ; Cook et al 1995 ; Jardí et al 1997 ; Denden et al 2009 ; Corda et al 2011	Not likely at increased risk ⁴	Corda et al 2011
PI*Q0 mattawa	1	Not likely at increased risk ¹	Cox & Levison 1988 ; Lara et al 2013	Not likely at increased risk ⁴	Cox & Levison 1988 ; Lara et al 2013
PI*Q0 clayton	1	Not likely at increased risk ¹	Rodriguez-Frias et al 2011	Not likely at increased risk ⁴	Clinical cases not reported
PI*M heerlen	1	Not likely at increased risk ¹	Kramps et al 1981 ; Poller et al 1999 ; Kalsheker et al 1992 ; Rodriguez et al 2002	Not likely at increased risk ⁴	Ferrarotti et al 2005 ; Balduyck et al 2014
PI*Z and PI*I	2	Slightly increased risk	Baur & Bencze 1987 ; Ferrarotti et al 2005 ; Corda et al 2006	Slightly increased risk ³	Baur & Bencze 1987 ; Mahadeva et al 1999 ; Ferrarotti et al 2005
PI*Z and PI*M procida	2	Unknown risk	Ferrarotti et al 2005 ; Lonardo et al 2002	Slightly increased risk ³	Ferrarotti et al 2005 ; Lonardo et al 2002
PI*Z and PI*M malton	2	Increased risk	Sproule et al 1983 ; Allen et al 1986 ; Ferrarotti et al 2005 ; Corda et al 2006 ; Suh-Lailam et al 2014 ; Joly et al 2015 ; Figueira Gonçalves et al 2017	Slightly increased risk ³	Sproule et al 1983 ; Ferrarotti et al 2005 ; Joly et al 2015 ; Figueira Gonçalves et al 2017

Genetic Result in the AlphalD™ At Home Genetic Health Risk Service Report	Number of variants	Reported Lung Disease Risk Category	References that support Lung Disease Risk Category	Reported Liver Disease Risk Category	References that support Liver Disease Risk Category
			Gonçalves et al 2017		
PI*Z and PI*S iiyama	2	Unknown risk	Clinical cases not reported	Slightly increased risk ³	Clinical cases not reported
PI*Z and PI*Q0 granite falls	2	Increased risk ²	Nukiwa et al 1987 ; Balduyck et al 2014 ; Pavičić et al 2019	Slightly increased risk ³	Clinical cases not reported
PI*Z and PI*Q0 west	2	Increased risk ²	Clinical cases not reported	Slightly increased risk ³	Clinical cases not reported
PI*Z and PI*Q0 bellingham	2	Increased risk ²	Poller et al 1990	Slightly increased risk ³	Poller et al 1990
PI*F and PI*Z	2	Slightly increased risk	Cockcroft et al 1981 ; Beckman et al 1984 ; Kelly et al 1989 ; Okayama et al 1991 ; Cook et al 1996 ; Sinden et al 2014 ; Franciosi et al 2019	Slightly increased risk ³	Kelly et al 1989 ; Sinden et al 2014 ; Franciosi et al 2019
PI*Z and PI*P lowell	2	Slightly increased risk	Holmes et al 1990 ; Ferrarotti et al 2005 ; Esteves-Brandão et al 2019 ; Bamforth & Kalsheker 1988	Slightly increased risk ³	Holmes et al 1990 ; Ferrarotti et al 2005 ; Bornhorst et al 2007 ; Bamforth & Kalsheker 1988
PI*Z and PI*Q0 mattawa	2	Increased risk ²	Balduyck et al 2014	Slightly increased risk ³	Clinical cases not reported
PI*Z and PI*Q0 clayton	2	Increased risk ²	Clinical cases not reported	Slightly increased risk ³	Clinical cases not reported
PI*Z and PI*M heerlen	2	Unknown risk	Klaassen et al 2001	Slightly increased risk ³	Clinical cases not reported
PI*S and PI*I	2	Unknown risk	Seri et al 1992 ; Huang et al 2017	Unknown risk	Clinical cases not reported
PI*S and PI*M procida	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*S and PI*M malton	2	Unknown risk	Figueira Gonçalves et al 2017	Slightly increased risk ³	Figueira Gonçalves et al 2017
PI*S and PI*S iiyama	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*S and PI*Q0 granite falls	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*S and PI*Q0 west	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*S and PI*Q0 bellingham	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*F and PI*S	2	Unknown risk	Cook et al 1996	Unknown risk	Clinical cases not reported
PI*S and PI*P lowell	2	Unknown risk	Jardí et al 1997 ; Balbi et al 2019	Unknown risk	Clinical cases not reported
PI*S and PI*Q0 mattawa	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported

Genetic Result in the AlphalD™ At Home Genetic Health Risk Service Report	Number of variants	Reported Lung Disease Risk Category	References that support Lung Disease Risk Category	Reported Liver Disease Risk Category	References that support Liver Disease Risk Category
PI*S and PI*Q0 clayton	2	Unknown risk	Rosenbaum et al 2017	Unknown risk	Rosenbaum et al 2017
PI*S and PI*M heerlen	2	Unknown risk	Kramps et al 1981	Unknown risk	Clinical cases not reported
PI*I and PI*I	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*I and PI*M procida	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*I and PI*M malton	2	Unknown risk	Clinical cases not reported	Slightly increased risk ³	Clinical cases not reported
PI*I and PI*S iiyama	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*I and PI*Q0 granite falls	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*I and PI*Q0 west	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*I and PI*Q0 bellingham	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*F and PI*I	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*I and PI*P lowell	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*I and PI*Q0 mattawa	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*I and PI*Q0 clayton	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*I and PI*M heerlen	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*F and PI*M procida	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*F and PI*M malton	2	Unknown risk	Clinical cases not reported	Slightly increased risk ³	Clinical cases not reported
PI*F and PI*S iiyama	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*F and PI*Q0 granite falls	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*F and PI*Q0 west	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*F and PI*Q0 bellingham	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*F and PI*F	2	Unknown risk	Sinden et al 2014	Not likely at increased risk ⁴	Sinden et al 2014
PI*F and PI*P lowell	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*F and PI*Q0 mattawa	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*F and PI*Q0 clayton	2	Unknown risk	Ringebach et al 2011	Not likely at increased risk ⁴	Ringebach et al 2011
PI*F and PI*M heerlen	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*M procida and PI*M procida	2	Unknown risk	Ferrarotti et al 2005	Not likely at increased risk ⁴	Ferrarotti et al 2005

Genetic Result in the AlphalD™ At Home Genetic Health Risk Service Report	Number of variants	Reported Lung Disease Risk Category	References that support Lung Disease Risk Category	Reported Liver Disease Risk Category	References that support Liver Disease Risk Category
PI*M malton and PI*M procida	2	Unknown risk	Clinical cases not reported	Slightly increased risk ³	Clinical cases not reported
PI*M procida and PI*S iiyama	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*M procida and PI*Q0 granite falls	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*M procida and PI*Q0 west	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*M procida and PI*Q0 bellingham	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*M procida and PI*P lowell	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*M procida and PI*Q0 mattawa	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*M procida and PI*Q0 clayton	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*M procida and PI*M heerlen	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*M malton and PI*M malton	2	Slightly increased risk	Curiel et al 1989b ; Ferrarotti et al 2005 ; Orrù et al 2005 ; Balduyck et al 2014 ; Joly et al 2015 ; Figueira Gonçalves et al 2017 ; Franciosi et al 2019 ; Reid et al 1987	Slightly increased risk ³	Curiel et al 1989b ; Janciauskienė et al 2004 ; Orrù et al 2005 ; Joly et al 2015 ; Figueira Gonçalves et al 2017 ; Reid et al 1987 ; Callea et al 2018
PI*M malton and PI*S iiyama	2	Unknown risk	Clinical cases not reported	Slightly increased risk ³	Clinical cases not reported
PI*M malton and PI*Q0 granite falls	2	Unknown risk	Clinical cases not reported	Slightly increased risk ³	Clinical cases not reported
PI*M malton and PI*Q0 west	2	Unknown risk	Clinical cases not reported	Slightly increased risk ³	Clinical cases not reported
PI*M malton and PI*Q0 bellingham	2	Unknown risk	Clinical cases not reported	Slightly increased risk ³	Clinical cases not reported
PI*M malton and PI*P lowell	2	Unknown risk	Clinical cases not reported	Slightly increased risk ³	Clinical cases not reported
PI*M malton and PI*Q0 mattawa	2	Unknown risk	Lara et al 2013 ; Balduyck et al 2014	Slightly increased risk ³	Lara et al 2013 ; Balduyck et al 2014
PI*M malton and PI*Q0 clayton	2	Unknown risk	Rodríguez-Frias et al 2011	Slightly increased risk ³	Clinical cases not reported
PI*M malton and PI*M heerlen	2	Unknown risk	Clinical cases not reported	Slightly increased risk ³	Clinical cases not reported
PI*S iiyama and PI*S iiyama	2	Slightly increased risk	Takabe et al 1992 ; Lomas et al 1993 ; Yuasa et al 1993	Unknown risk	Takabe et al 1992 ; Yuasa et al 1993
PI*Q0 granite falls and PI*S iiyama	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*S iiyama and PI*Q0 west	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*S iiyama and PI*Q0 bellingham	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported

Genetic Result in the AlphalD™ At Home Genetic Health Risk Service Report	Number of variants	Reported Lung Disease Risk Category	References that support Lung Disease Risk Category	Reported Liver Disease Risk Category	References that support Liver Disease Risk Category
PI*P lowell and PI*S iiyama	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*S iiyama and PI*Q0 mattawa	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*S iiyama and PI*Q0 clayton	2	Unknown risk	Ko et al 2011; Miyahara et al 2001	Unknown risk	Ko et al 2011; Miyahara et al 2001
PI*M heerlen and PI*S iiyama	2	Unknown risk	Clinical cases not reported	Unknown risk	Clinical cases not reported
PI*Q0 granite falls and PI*Q0 granite falls	2	Increased risk ²	Holmes et al 1989; Hubbard et al 1989; Balbi et al 2019	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 granite falls and PI*Q0 west	2	Increased risk ²	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 bellingham and PI*Q0 granite falls	2	Increased risk ²	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*P lowell and PI*Q0 granite falls	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 granite falls and PI*Q0 mattawa	2	Increased risk ²	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 clayton and PI*Q0 granite falls	2	Increased risk ²	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*M heerlen and PI*Q0 granite falls	2	Unknown risk	Poller et al 1999	Not likely at increased risk ⁴	Poller et al 1999
PI*Q0 west and PI*Q0 west	2	Increased risk ²	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 bellingham and PI*Q0 west	2	Increased risk ²	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*P lowell and PI*Q0 west	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 mattawa and PI*Q0 west	2	Increased risk ²	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 clayton and PI*Q0 west	2	Increased risk ²	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*M heerlen and PI*Q0 west	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 bellingham and PI*Q0 bellingham	2	Increased risk ²	Satoh et al 1988; Cook et al 1994; Garver et al 1986; Fregonese et al 2008	Not likely at increased risk ⁴	Cook et al 1994; Garver et al 1986
PI*P lowell and PI*Q0 bellingham	2	Unknown risk	Cook et al 1994; Cook et al 1995	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 bellingham and PI*Q0 mattawa	2	Increased risk ²	Curiel et al 1989a	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 bellingham and PI*Q0 clayton	2	Increased risk ²	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*M heerlen and PI*Q0 bellingham	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*P lowell and PI*P lowell	2	Unknown risk	Faber et al 1989	Not likely at increased risk ⁴	Clinical cases not reported

Genetic Result in the AlphaID™ At Home Genetic Health Risk Service Report	Number of variants	Reported Lung Disease Risk Category	References that support Lung Disease Risk Category	Reported Liver Disease Risk Category	References that support Liver Disease Risk Category
PI*P lowell and PI*Q0 mattawa	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*P lowell and PI*Q0 clayton	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*P lowell and PI*M heerlen	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 mattawa and PI*Q0 mattawa	2	Increased risk ²	Cox & Levinson 1988	Not likely at increased risk ⁴	Cox & Levinson 1988
PI*Q0 clayton and PI*Q0 mattawa	2	Increased risk ²	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*M heerlen and PI*Q0 mattawa	2	Unknown risk	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*Q0 clayton and PI*Q0 clayton	2	Increased risk ²	Clinical cases not reported	Not likely at increased risk ⁴	Clinical cases not reported
PI*M heerlen and PI*Q0 clayton	2	Unknown risk	Brantly et al 1997	Not likely at increased risk ⁴	Clinical cases not reported
PI*M heerlen and PI*M heerlen	2	Slightly increased risk	Kramps et al 1981; Fregonese et al 2008	Not likely at increased risk ⁴	Kramps et al 1981
Variant not determined	N/A	Not determined risk	N/A	Not determined risk	N/A

¹All carriers (except ever-smokers PI*M/PI*Z) are reported as “Not likely at increased risk” for lung disease risk even if insufficient number of clinical cases have been reported.

²PI*Z/PI*Null and PI*Null/PI*Null combinations are assigned as “Increased risk” for lung disease risk even if insufficient number of clinical cases have been reported. Null variants: PI*Q0 granite falls, PI*Q0 west, PI*Q0 bellingham, PI*Q0 mattawa and PI*Q0 clayton.

³Individuals with at least one PI*Z or PI*M malton severe accumulation variants and their combinations are reported as “Slightly Increased risk” for liver disease risk even if insufficient number of clinical cases have been reported.

⁴All the combinations including two non-accumulation alleles are reported as “Not likely at Increased risk” for liver disease risk even if insufficient number of clinical cases have been reported. Non-accumulation alleles: PI*M, PI*M procida, PI*Q0 granite falls, PI*Q0 west, PI*Q0 bellingham, PI*F, PI*P lowell, PI*Q0 mattawa, PI*Q0 clayton, and PI*M heerlen.

- **Abboud et al 1979:** Abboud RT, Rushton JM, Grzybowski S. Interrelationships between neutrophil elastase, serum alpha, -antitrypsin, lung function and chest radiography in patients with chronic airflow obstruction. Am Rev Respir Dis. 1979 Jul;120(1):31-40. doi: 10.1164/arrd.1979.120.1.31. PMID: 313728.
- **Allen et al 1986:** Allen MB, Ward AM, Perks WH. Alpha 1 antitrypsin deficiency due to MMaltonZ phenotype: case report and family study. Thorax. 1986 Jul;41(7):568-70. PubMed PMID: 3491442; PubMed Central PMCID: PMC460392.
- **ATS/ERS 2003:** American Thoracic Society/European Respiratory Society statement: standards for the diagnosis and management of individuals with alpha-1 antitrypsin deficiency. American Thoracic Society; European Respiratory Society. Am J Respir Crit Care Med. 2003 Oct 1;168(7):818-900.
- **Arnaud et al 1977:** Arnaud P, Chapuis-Cellier C, Vittoz P, Creyssel R. Alpha-1-antitrypsin phenotypes in Lyon, France. Hum Genet. 1977;39(1):63-68. doi:10.1007/BF00273153

- **Arnaud et al 1978:** Arnaud P, Fudenberg H.H, Chapuis-Cellier C, Vittoz, P. Genetic polymorphism of serum alpha-1-protease inhibitor (alpha-1-antitrypsin): Pi I, a deficient allele of the Pi system. In: *The Journal of Laboratory and Clinical Medicine*. (The Journal of Laboratory and Clinical Medicine, August 1978, 92(2):177-184).
- **Balbi et al 2019:** Balbi B, Sangiorgi C, Gnementi I, Ferrarotti I, Vallese D, Paracchini E, Delle Donne L, Corda L, Baderna P, Corsico A, Carone M, Brun P, Cappello F, Ricciardolo FL, Ruggeri P, Mumby S, Adcock IM, Caramori G, Di Stefano A. Bacterial load and inflammatory response in sputum of alpha-1 antitrypsin deficiency patients with COPD. *Int J Chron Obstruct Pulmon Dis.* 2019 Aug 21;14:1879-1893. doi: 10.2147/COPD.S207203.
- **Baldyuck et al 2014:** Baldyuck M, Odou MF, Zerimech F, Porchet N, Lafitte JJ, Maitre B. Diagnosis of alpha-1 antitrypsin deficiency: modalities, indications and diagnosis strategy. *Rev Mal Respir.* 2014 Oct;31(8):729-45. doi: 10.1016/j.rmr.2014.06.001. Review. PubMed PMID: 25391508.
- **Bamforth & Kalsheker 1988:** Bamforth F, Kalsheker N. Alpha1 antitrypsin deficiency due to Pi null: clinical presentation and evidence for molecular heterogeneity. *J Med Genet.* 1988. 25, 83-87.
- **Bartmann et al 1985:** Bartmann K, Fooke-Achterrath M, Koch G, Nagy I, Schütz I, Weis E, Zierski M. Heterozygosity in the Pi-system as a pathogenetic cofactor in chronic obstructive pulmonary disease (COPD). *Eur J Respir Dis.* 1985 Apr;66(4):284-96. PMID: 3874784.
- **Baur & Bencze 1987:** Baur X, Bencze K. Study of familial alpha-1-proteinase inhibitor deficiency including a rare proteinase inhibitor phenotype (IZ). I. Alpha-1-phenotyping and clinical investigations. *Respiration.* 1987;51(3):188-95. PubMed PMID: 3496639.
- **Beckman et al 1984:** Beckman G, Stjernberg NL, Eklund A. Is the PiF allele of alpha 1-antitrypsin associated with pulmonary disease? *Clin Genet.* 1984 25(6):491-5. PubMed PMID: 6610506.
- **Bell et al 1990:** Bell H, Schrumpf E, Fagerhol MK. Heterozygous MZ alpha-1-antitrypsin deficiency in adults with chronic liver disease. *Scand J Gastroenterol.* 1990;25(8):788-792. doi:10.3109/00365529008999216
- **Blenkinsopp & Haffenden 1977:** Blenkinsopp WK, Haffenden GP. Alpha-1-antitrypsin bodies in the liver. *J Clin Pathol.* 1977;30(2):132-137. doi:10.1136/jcp.30.2.132
- **Bornhorst et al 2007:** Bornhorst JA, Calderon FR, Procter M, Tang W, Ashwood ER, Mao R. Genotypes and serum concentrations of human alpha-1-antitrypsin "P" protein variants in a clinical population. *J Clin Pathol.* 2007 60(10):1124-8. PubMed PMID: 17906067.
- **Bornhorst et al 2013:** Bornhorst JA, Greene DN, Ashwood ER, Grenache DG. a1-Antitrypsin phenotypes and associated serum protein concentrations in a large clinical population. *Chest.* 2013 Apr;143(4):1000-1008. doi: 10.1378/chest.12-0564. PubMed PMID: 23632999.
- **Brantly et al 1997:** Brantly M, Lee JH, Hildesheim J, Uhm CS, Prakash UB, Staats BA, Crystal RG, Hildesheim J. alpha1-antitrypsin gene mutation hot spot associated with the formation of a retained and degraded null variant. *Am J Respir Cell Mol Biol.* 1997. Mar;16(3):225-31.
- **Callea et al 2018:** Callea F, Giovannoni I, Francalanci P, Boldrini R, Faa G, Medicina D, Nobili V, Desmet VJ, Ishak K, Seyama K, Bellacchio E. Mineralization of alpha-1-antitrypsin inclusion bodies in Mmalton alpha-1-antitrypsin deficiency. *Orphanet J Rare Dis.* 2018 May 16;13(1):79. doi: 10.1186/s13023-018-0821-7.
- **Canva et al 2001:** Canva V, Piotte S, Aubert JP, Porchet N, Lecomte-Houcke M, Huet G, Zenjari T, Roumilhac D, Pruvot FR, Degand P, Paris JC, Balduyck M. Heterozygous

- M3Mmalton alpha1-antitrypsin deficiency associated with end-stage liver disease: case report and review. Clin Chem. 2001 Aug;47(8):1490-6. PubMed PMID: 11468249.
- **Carlson et al 1981:** Carlson J, Eriksson S, Hägerstrand I. Intra- and extracellular alpha 1-antitrypsin in liver disease with special reference to Pi phenotype. J Clin Pathol. 1981;34(9):1020-1025. doi:10.1136/jcp.34.9.1020
 - **Carroll et al 2011:** Carroll TP, O'Connor CA, Floyd O, McPartlin J, Kelleher DP, O'Brien G, Dimitrov BD, Morris VB, Taggart CC, McElvaney NG. The prevalence of alpha-1 antitrypsin deficiency in Ireland. Respir Res. 2011 Jul 13;12:91. doi: 10.1186/1465-9921-12-91. PubMed PMID: 21752289.
 - **Chan-Yeung et al 1978:** Chan-Yeung M, Ashley MJ, Corey P, Maledy H. Pi phenotypes and the prevalence of chest symptoms and lung function abnormalities in workers employed in dusty industries. Am Rev Respir Dis. 1978 Feb;117(2):239-45. doi: 10.1164/arrd.1978.117.2.239. PMID: 305739.
 - **Cockcroft et al 1981:** Cockcroft DW, Tennent RK, Horne SL. Pulmonary emphysema associated with the FZ alpha 1-antitrypsin phenotype. Can Med Assoc J. 1981 15;124(6):737-42. PubMed PMID: 6970613.
 - **Cook et al 1994:** Cook L, Janus ED, Brenton S, Tai E, Burdon J. Absence of alpha-1-antitrypsin (Pi Null Bellingham) and the early onset of emphysema. Aust N Z J Med. 1994 Jun;24(3):263-9. PubMed PMID: 7980208.
 - **Cook et al 1995:** Cook L, Burdon J, Brenton S, Janus ED, Knight K. Alpha-1-antitrypsin PLowell: a normally functioning variant present in low concentration. Aust N Z J Med. 1995 Dec;25(6):695-7. PubMed PMID: 8770333.
 - **Cook et al 1996:** Cook L, Burdon JG, Brenton S, Knight KR, Janus ED. Kinetic characterisation of alpha-1-antitrypsin F as an inhibitor of human neutrophil elastase. Pathology. 1996 28(3):242-7. PubMed PMID: 8912354.
 - **Corda et al 2006:** Corda L, Bertella E, Pini L, Pezzini A, Medicina D, Boni E, Guerini M, Trivella S, Grassi V, Tantucci C. Diagnostic flow chart for targeted detection of alpha1-antitrypsin deficiency. Respir Med. 2006 Mar;100(3):463-70. PubMed PMID: 16043335
 - **Corda et al 2011:** Corda L, Medicina D, La Piana GE, Bertella E, Moretti G, Bianchi L, Pinelli V, Savoldi G, Baiardi P, Facchetti F, Gatta N, Annesi-Maesano I, Balbi B. Population genetic screening for alpha1-antitrypsin deficiency in a high-prevalence area. Respiration. 2011;82(5):418-25. doi: 10.1159/000325067. PubMed PMID: 21474916.
 - **Cox & Levison 1988:** Cox DW, Levison H. Emphysema of early onset associated with a complete deficiency of alpha-1-antitrypsin (null homozygotes). Am Rev Respir Dis. 1988 137(2):371-5. PubMed PMID: 3257661.
 - **Curiel et al 1989a:** Curiel D, Brantly M, Curiel E, Stier L, Crystal RG. Alpha 1-antitrypsin deficiency caused by the alpha 1-antitrypsin Null mattawa gene. An insertion mutation rendering the alpha 1-antitrypsin gene incapable of producing alpha 1-antitrypsin. J Clin Invest. 1989. Apr; 83(4):1144-52.
 - **Curiel et al 1989b:** Curiel DT, Holmes MD, Okayama H, Brantly ML, Vogelmeier C, Travis WD, Stier LE, Perks WH, Crystal RG. Molecular basis of the liver and lung disease associated with the alpha 1-antitrypsin deficiency allele Mmalton. J Biol Chem. 1989. 15; 264(23):13938-45. (Curiel 1989b)
 - **Dahl et al 2002:** Dahl M, Tybjaerg-Hansen A, Lange P, Vestbo J, Nordestgaard BG. Change in lung function and morbidity from chronic obstructive pulmonary disease in alpha1-antitrypsin MZ heterozygotes: A longitudinal study of the general population. Ann Intern Med. 2002 Feb 19;136(4):270-9. doi: 10.7326/0003-4819-136-4-200202190-00006. PMID: 11848724.
 - **de Serres & Blanco 2014:** de Serres F, Blanco I. Role of alpha-1 antitrypsin in human health and disease. J Intern Med. 2014 Oct;276(4):311-35. doi: 10.1111/joim.12239.

- **DeMeo & Silverman 2004:** DeMeo DL, Silverman EK. Alpha1-antitrypsin deficiency. 2: genetic aspects of alpha(1)-antitrypsin deficiency: phenotypes and genetic modifiers of emphysema risk. *Thorax*. 2004 Mar;59(3):259-64. Review. PubMed PMID: 14985567
- **Denden et al 2009:** Denden S, Zorzetto M, Amri F, Knani J, Ottaviani S, Scabini R, Gorrini M, Ferrarotti I, Campo I, Chibani JB, Khelil AH, Luisetti M. Screening for Alpha 1 antitrypsin deficiency in Tunisian subjects with obstructive lung disease: a feasibility report. *Orphanet J Rare Dis*. 2009 Apr 15;4:12. doi: 10.1186/1750-1172-4-12. PubMed PMID: 19368725.
- **Duk et al 2016:** Duk K, Zdral A, Szumna B, Roży A, Chorostowska-Wynimko J. Frequency of Rare Alpha-1 Antitrypsin Variants in Polish Patients with Chronic Respiratory Disorders. *Adv Exp Med Biol*. 2016;910:47-53. doi: 10.1007/5584_2016_213. PubMed PMID: 26987331
- **Eigenbrodt et al 1997:** Eigenbrodt ML, McCashland TM, Dy RM, Clark J, Galati J. Heterozygous alpha 1-antitrypsin phenotypes in patients with end stage liver disease. *Am J Gastroenterol*. 1997;92(4):602-607.
- **Elzouki & Eriksson 1996:** Elzouki AN, Eriksson S. Risk of hepatobiliary disease in adults with severe alpha 1-antitrypsin deficiency (PiZZ): is chronic viral hepatitis B or C an additional risk factor for cirrhosis and hepatocellular carcinoma?. *Eur J Gastroenterol Hepatol*. 1996;8(10):989-994.
- **Eriksson et al 1986:** Eriksson S, Carlson J, Velez R. Risk of cirrhosis and primary liver cancer in alpha 1-antitrypsin deficiency. *N Engl J Med*. 1986 Mar 20;314(12):736-9. doi: 10.1056/NEJM198603203141202. PMID: 3485248.
- **Esteves Brandão et al 2019:** Esteves Brandão M, Conde B, Seixas S, Clotilde Silva J, Fernandes A. Pulmonary Emphysema in a Child With Alpha-1 Antitrypsin Deficiency: Evaluation of 2 Years of Intravenous Augmentation Therapy. *Arch Bronconeumol*. 2019 Sep;55(9):502-504.
- **Faber et al 1989:** Faber JP, Weidinger S, Goedde HW, Ole K. The deficient alpha-1-antitrypsin phenotype PI P is associated with an A-to-T transversion in exon III of the gene. *Am J Hum Genet*. 1989 Jul;45(1):161-3.
- **Fagerhol & Hauge 1969:** Fagerhol MK, Hauge HE. Serum Pi types in patients with pulmonary diseases. *Acta Allergol*. 1969 May;24(2):107-14. doi: 10.1111/j.1398-9995.1969.tb03760.x. PMID: 4191350.
- **Ferrarotti et al 2005:** Ferrarotti I, Baccheschi J, Zorzetto M, Tinelli C, Corda L, Balbi B, Campo I, Pozzi E, Faa G, Coni P, Massi G, Stella G, Luisetti M. Prevalence and phenotype of subjects carrying rare variants in the Italian registry for alpha1-antitrypsin deficiency. *J Med Genet*. 2005 Mar;42(3):282-7. PubMed PMID: 15744045.
- **Figueira Gonçalves et al 2017:** Figueira Gonçalves JM, Martínez Bugallo F, Díaz Pérez D, Martín Martínez MD, García-Talavera I, Pitti Pérez R. Clinical manifestations of the Mmalton alpha-1 antitrypsin deficiency variant. *Pulmonology*. 2017;24:48-49.
- **Franciosi et al 2019:** Franciosi AN, Carroll TP, McElvaney NG. Pitfalls and caveats in α1-antitrypsin deficiency testing: a guide for clinicians. *Lancet Respir Med*. 2019 Dec;7(12):1059-1067.
- **Fregonese et al 2008:** Fregonese L, Stolk J, Frants RR, Veldhuisen B. Alpha-1 antitrypsin Null mutations and severity of emphysema. *Respir Med*. 2008 Jun;102(6):876-84. doi: 10.1016/j.rmed.2008.01.009. PubMed PMID: 18353624.
- **Garver et al 1986:** Garver RI Jr, Mornex JF, Nukiwa T, Brantly M, Courtney M, LeCocq JP, Crystal RG. Alpha 1-antitrypsin deficiency and emphysema caused by homozygous inheritance of non-expressing alpha 1-antitrypsin genes. *N Engl J Med*. 1986 Mar 20;314(12):762-6.

- **Graham et al 2015:** Graham RP, Dina MA, Howe SC, Butz ML, Willkomm KS, Murray DL, Snyder MR, Rumilla KM, Halling KC, Highsmith WE Jr. SERPINA1 Full-Gene Sequencing Identifies Rare Mutations Not Detected in Targeted Mutation Analysis. *J Mol Diagn.* 2015 Nov;17(6):689-94. doi: 10.1016/j.jmoldx.2015.07.002. PubMed PMID: 26321041.
- **Graziadei et al 1998:** Graziadei IW, Joseph JJ, Wiesner RH, Therneau TM, Batts KP, Porayko MK. Increased risk of chronic liver failure in adults with heterozygous alpha1-antitrypsin deficiency. *Hepatology.* 1998;28(4):1058-1063. doi:10.1002/hep.510280421
- **Gulsvik et al 1979:** Gulsvik A, Fagerhol MK. Alpha 1-antitrypsin phenotypes and obstructive lung disease in the city of Oslo. *Scand J Respir Dis.* 1979 Oct;60(5):267-74. PMID: 316573.
- **Hamesch et al 2019:** Hamesch K, Mandorfer M, Pereira VM, Moeller LS, Pons M, Dolman GE, Reichert MC, Schneider CV, Woditsch V, Voss J, Lindhauer C, Fromme M, Spivak I, Guldiken N, Zhou B, Arslanow A, Schaefer B, Zoller H, Aigner E, Reiberger T, Wetzel M, Siegmund B, Simões C, Gaspar R, Maia L, Costa D, Bento-Miranda M, van Helden J, Yagmur E, Bzdok D, Stolk J, Gleiber W, Knipel V, Windisch W, Mahadeva R, Bals R, Koczulla R, Barrecheguren M, Miravitles M, Janciauskiene S, Stickel F, Lammert F, Liberal R, Genesca J, Griffiths WJ, Trauner M, Krag A, Trautwein C, Strnad P; European Alpha1-Liver Study Group.. Liver Fibrosis and Metabolic Alterations in Adults With alpha-1-antitrypsin Deficiency Caused by the Pi*ZZ Mutation. *Gastroenterology.* 2019;157(3):705-719.e18. doi:10.1053/j.gastro.2019.05.013
- **Holmes et al 1989:** Holmes M, Curiel D, Brantly M, Crystal RG. Characterization of the intracellular mechanism causing the alpha-1-antitrypsin Nullgranite falls deficiency state. *Am Rev Respir Dis.* 1989 Dec;140(6):1662-7.
- **Holmes et al 1990:** Holmes MD, Brantly ML, Crystal RG. Molecular analysis of the heterogeneity among the P-family of alpha-1-antitrypsin alleles. *Am Rev Respir Dis.* 1990 142(5):1185-92. PubMed PMID: 2240842.
- **Huang et al 2017:** Huang AC, Perez C, Felipez L, Chamyan G, Finegold MJ, Hernandez E. A Rare Phenotype of Alpha-1-Antitrypsin Deficiency Owing to PI*IS in a Newborn With Liver Disease. *J Pediatr Gastroenterol Nutr.* 2017 Nov;65(5):e112-e113. doi: 10.1097. PubMed PMID: 28837509.
- **Hubbard et al 1989:** Hubbard RS, McElvaney NG, Sellers SE, Healy JT, Czerski DB, Crystal RG. Recombinant DNA-produced alpha 1-antitrypsin administered by aerosol augments lower respiratory tract antineutrophil elastase defenses in individuals with alpha 1-antitrypsin deficiency. *J Clin Invest.* 1989 Oct; 84(4): 1349–1354.
- **Janciauskiene et al 2004:** Janciauskiene S, Eriksson S, Callea F, Mallya M, Zhou A, Seyama K, Hata S, Lomas DA. Differential detection of PAS-positive inclusions formed by the Z, Siyama, and Mmalton variants of alpha1-antitrypsin. *Differential Hepatology.* 2004 Nov;40(5):1203-10.
- **Jardí et al 1997:** Jardí R, Rodríguez-Frías F, Casas F, Cotrina M, Vidal R, Miravitles M, Pascual C. [Molecular characterization of two variants of alpha-1-antitrypsin deficiency: PI Mpalermo and PI Plovel]. *Med Clin (Barc).* 1997 11;109(12):463-6. Spanish. PubMed PMID: 9441182.
- **Joly et al 2015:** Joly P, Guillaud O, Hervieu V, Francina A, Mornex JF, Chapuis-Cellier C. Clinical heterogeneity and potential high pathogenicity of the Mmalton Alpha 1 antitrypsin allele at the homozygous, compound heterozygous and heterozygous states. *Orphanet J Rare Dis.* 2015 Oct 7;10:130. doi: 10.1186/s13023-015-0350-6.
- **Kalsheker et al 1992:** Kalsheker N, Hayes K, Weidinger S, Graham A. What is Pi (proteinase inhibitor) null or PiQ0?: a problem highlighted by the alpha 1 antitrypsin Mheerlen mutation. *J Med Genet.* 1992 Jan;29(1):27-9. PubMed PMID: 1552539; PubMed Central PMCID: PMC1015817.

- **Kelly et al 1989:** Kelly CP, Tyrrell DN, McDonald GS, Whitehouse DB, Prichard JS. Heterozygous FZ alpha 1 antitrypsin deficiency associated with severe emphysema and hepatic disease: case report and family study. Thorax. 1989 Sep;44(9):758-9. PubMed PMID: 2588214.
- **Klaassen et al 2001:** Klaassen CH, de Metz M, van Aarsen Y, Janssen J. alpha(1)-Antitrypsin deficiency as a result of compound heterozygosity for the Z and M(Heerlen) alleles. Clin Chem. 2001 May;47(5):978-9.
- **Klayton et al 1975:** Klayton R, Fallat R, Cohen AB. Determinants of chronic obstructive pulmonary disease in patients with intermediate levels of alpha-antitrypsin. Am Rev Respir Dis. 1975 Jul;112(1):71-5. doi: 10.1164/arrd.1975.112.1.71. PMID: 238443.
- **Ko et al 2011:** Ko DH, Chang HE, Song SH, Yoon H, Park KU, Song J. Identification of compound heterozygous mutation in a Korean patient with alpha 1-antitrypsin deficiency. Korean J Lab Med. 2011 Oct;31(4):294-7. doi: 10.3343/kjlm.2011.31.4.294. PubMed PMID: 22016686; PubMed Central PMCID: PMC3190011.
- **Kramps et al 1981:** Kramps JA, Brouwers JW, Maesen F, Dijkman JH. PiMheerlen, alpha PiM allele resulting in very low alpha 1-antitrypsin serum levels. Hum Genet. 1981;59(2):104-7. PubMed PMID: 6976926.
- **Kueppers et al 1977:** Kueppers F, Miller RD, Gordon H, Hepper NG, Offord K. Familial prevalence of chronic obstructive pulmonary disease in a matched pair study. Am J Med. 1977 Sep;63(3):336-42. doi: 10.1016/0002-9343(77)90270-4. PMID: 302643.
- **Kwok et al 2004:** Kwok JS, Lawton JW, Yew WW, Chau CH, Lee J, Wong PC. Protease inhibitor phenotypes and serum alpha-1-antitrypsin levels in patients with COPD: a study from Hong Kong. Respirology. 2004 Jun;9(2):265-70. PubMed PMID: 15182280.
- **Lang et al 2005:** Lang T, Mühlbauer M, Strobelt M, Weidinger S, Hadorn HB. Alpha-1-antitrypsin deficiency in children: liver disease is not reflected by low serum levels of alpha-1-antitrypsin - a study on 48 pediatric patients. Eur J Med Res. 2005 Dec 7;10(12):509-14. PubMed PMID: 16356865.
- **Lara et al 2013:** Lara B, Martínez-Delgado B, Torres ML, Marín-Arguedas S, Bustamante A, Miravitlles M. Alpha-1-antitrypsin deficiency associated with the Mattawa variant. Arch Bronconeumol. 2013. Dec;49(12):548-50.
- **Laubach et al 1993:** Laubach VE, Ryan WJ, Brantly M. Characterization of a human alpha 1-antitrypsin null allele involving aberrant mRNA splicing. Hum Mol Genet. 1993 Jul;2(7):1001-5.
- **Lieberman et al 1986:** Lieberman J, Winter B, Sastre A. Alpha 1-antitrypsin Pi-types in 965 COPD patients. Chest. 1986 Mar;89(3):370-3. doi: 10.1378/chest.89.3.370. PMID: 3485034.
- **Lochon et al 1978:** Lochon B, Vercaigne D, Lochon C, Fournier M, Martin JP, Derenne JP, Pariente R. Emphysème pan-lobulaire: relations avec le taux d'alpha 1-antitrypsine sérique, le phénotype Pi et le système H.L.A [Pan-lobular emphysema: relationship with serum alpha-1-antitrypsin levels, Pi phenotype and the HLA system (author's transl)]. Nouv Presse Med. 1978 Apr 8;7(14):1167-70. French. PMID: 307224.
- **Lomas et al 1993:** Lomas DA, Finch JT, Seyama K, Nukiwa T, Carrell RW. Alpha 1-antitrypsin Siiyama (Ser53-->Phe). Further evidence for intracellular loop-sheet polymerization. J Biol Chem. 1993 Jul 25;268(21):15333-5. PubMed PMID: 8340361.
- **Lonardo et al 2002:** Lonardo A, Medicina D, Leonelli M, Bagni A, Callea F. Intestinal Wegener's granulomatosis in a patient with severe alpha-1-antitrypsin deficiency resulting from a unique combination of two deficiency alleles (PiZ and PiMProcida). Eur J Gastroenterol Hepatol. 2002 Dec;14(12):1389-92.
- **Mahadeva et al 1999:** Mahadeva R, Chang WS, Dafforn TR, Oakley DJ, Foreman RC,

- Calvin J, Wight DG, Lomas DA. Heteropolymerization of S, I, and Z alpha1-antitrypsin and liver cirrhosis. *J Clin Invest.* 1999 Apr;103(7):999-1006. PubMed PMID: 10194472
- **Matzen et al 1977:** Matzen RN, Bader PI, Block WD. alpha1-Antitrypsin deficiency in clinic patients. *Ann Clin Res.* 1977 Apr;9(2):88-92. PMID: 302107.
 - **McElvaney et al 2020:** McElvaney GN, Sandhaus RA, Miravitles M, Turino GM, Seersholtz N, Wencker M, Stockley RA. Clinical considerations in individuals with Alpha-1 Antitrypsin PI*SZ genotype. *Eur Respir J.* 2020 Mar 12. doi: 10.1183/13993003.02410-2019.
 - **Meira et al 2018:** Meira L, Boaventura R, Seixas S, Sucena M. Alpha-1 Antitrypsin Deficiency Detection in a Portuguese Population. *COPD.* 2018 15(1):4-9. doi: 10.1080/15412555.2017.1414779. PubMed PMID: 29393705.
 - **Miyahara et al 2001:** Miyahara N, Seyama K, Sato T, Fukuchi Y, Eda R, Takeyama H, Harada M. Compound heterozygosity for alpha-1-antitrypsin (S(iiyama) and Q0(clayton)) in an Oriental patient. *Intern Med.* 2001 Apr;40(4):336-40. PubMed PMID: 11334395.
 - **Molloy et al 2014:** Molloy K, Hersh CP, Morris VB, Carroll TP, O'Connor CA, Lasky-Su JA, Greene CM, O'Neill SJ, Silverman EK, McElvaney NG. Clarification of the risk of chronic obstructive pulmonary disease in α 1-antitrypsin deficiency PiMZ heterozygotes. *Am J Respir Crit Care Med.* 2014 Feb 15;189(4):419-27. doi: 10.1164/rccm.201311-1984OC.
 - **Montealegre et al 2006:** Montealegre F, Delgado A, Toro A, Vargas W, Chardon D, Bayona M, Campbell E. Alfa 1 antitrypsin and protease levels in Puerto Rican asthmatics: a pilot study. *P R Health Sci J.* 2006 Jun;25(2):117-25. PubMed PMID: 17203708.
 - **Nukiwa et al 1987:** Nukiwa T, Takahashi H, Brantly M, Courtney M, Crystal RG. Alpha 1-Antitrypsin nullGranite Falls, a nonexpressing alpha 1-antitrypsin gene associated with a frameshift to stop mutation in a coding exon. *J Biol Chem.* 1987. Sep 5; 262(25):11999-2004.
 - **Okayama et al 1991:** Okayama H, Brantly M, Holmes M, Crystal RG. Characterization of the molecular basis of the alpha 1-antitrypsin F allele. *Am J Hum Genet.* 1991. Jun; 48(6):1154-8.
 - **Orrù et al 2005:** Orrù G, Faa G, Pillai S, Pilloni L, Montaldo C, Pusceddu G, Piras V, Coni P. Rapid PCR real-time genotyping of M-Malton alpha1-antitrypsin deficiency alleles by molecular beacons. *Diagn Mol Pathol.* 2005 Dec;14(4):237-42. PubMed PMID: 16319694.
 - **Ostrow et al 1978:** Ostrow DN, Manfreda J, Dorman T, Cherniack RM. Alpha1-antitrypsin phenotypes and lung function in a moderately polluted northern Ontario community. *Can Med Assoc J.* 1978 Mar 18;118(6):669-72. PMID: 306869; PMCID: PMC1818022.
 - **Pavičić et al 2019:** Pavičić T, Ćelap I, Njegovan M, Tešija Kuna A, Štefanović M. α -1 Antitrypsin Genotype-Phenotype Discrepancy in a 42-Year-Old Man Who Carries the Null-Allele. *Lab Med.* 2019 Oct 4. pii: lmz059. doi: 10.1093/labmed/lmz059.
 - **Poller et al 1990:** Poller W, Faber JP, Olek K. Highly variable clinical course in severe alpha 1-antitrypsin deficiency--use of polymerase chain reaction for the detection of rare deficiency alleles. *Klin Wochenschr.* 1990 Sep 3;68(17):857-63.
 - **Poller et al 1999:** Poller W, Merklein F, Schneider-Rasp S, Haack A, Fechner H, Wang H, Anagnostopoulos I, Weidinger S. Molecular characterisation of the defective alpha 1-antitrypsin alleles PI PI*Mwurzburg (Pro369Ser), PI*M heerlen (Pro369Leu), and PI*Q0 Lisbon (Thr68Ile). *Eur J Hum Genet.* 1999. Apr;7(3):321-31.
 - **Prins et al 2008:** Prins J, van der Meijden BB, Kraaijenhagen RJ, Wielders JP. Inherited chronic obstructive pulmonary disease: new selective-sequencing workup for alpha1-antitrypsin deficiency identifies 2 previously unidentified null alleles. *Clin Chem.* 2008 54(1):101-7. PubMed PMID: 18024524.

- **Propst et al 1992:** Propst T, Propst A, Dietze O, Judmaier G, Braunsteiner H, Vogel W. High prevalence of viral infection in adults with homozygous and heterozygous alpha 1-antitrypsin deficiency and chronic liver disease. *Ann Intern Med.* 1992;117(8):641-645. doi:10.7326/0003-4819-117-8-641
- **Reid et al 1987:** Reid CL, Wiener GJ, Cox DW, Richter JE, Geisinger KR. Diffuse hepatocellular dysplasia and carcinoma associated with the Mmalton variant of alpha 1-antitrypsin. *Gastroenterology.* 1987;93(1):181–7.
- **Ringenbach et al 2011:** Ringenbach MR, Banta E, Snyder MR, Craig TJ, Ishmael FT. A challenging diagnosis of alpha-1-antitrypsin deficiency: identification of a patient with a novel F/Null phenotype. *Allergy Asthma Clin Immunol.* 2011. Nov 13;7(1):18.
- **Rodriguez et al 2002:** Rodriguez F, Jardí R, Costa X, Cotrina M, Galimany R, Vidal R, Miravitles M. Rapid screening for alpha1-antitrypsin deficiency in patients with chronic obstructive pulmonary disease using dried blood specimens. *Am J Respir Crit Care Med.* 2002 Sep 15;166(6):814-7.
- **Rodríguez-Frias et al 2011:** Rodríguez-Frias F, Vila B, Homs M, Vidal R, Calpe JL, Jardi R. Diagnosis of Alpha-1 Antitrypsin Deficiency: Limitations of Rapid Diagnostic Laboratory Tests. *Arch Bronconeumol.* 2011. Aug;47(8):415-7.
- **Rodríguez-Frias et al 2012:** Rodriguez-Frias F, Miravitles M, Vidal R, Camos S, Jardi R. Rare alpha-1-antitrypsin variants: are they really so rare? *Ther Adv Respir Dis.* 2012 Apr;6(2):79-85. doi: 10.1177/1753465811434320. PubMed PMID: 22291048.
- **Rosenbaum et al 2017:** Rosenbaum EM, Chapaton-Rivard E, Overdorf C. Alpha 1 Antitrypsin Deficiency, Two Cases of Heterozygous S and Clayton Null Alleles. Volume 2 Number 2 Winter, 2017 Pages 22-28
- **Sandford et al 1999:** Sandford AJ, Weir TD, Spinelli JJ, Paré PD. Z and S mutations of the alpha1-antitrypsin gene and the risk of chronic obstructive pulmonary disease. *Am J Respir Cell Mol Biol.* 1999 Feb;20(2):287-91. doi: 10.1165/ajrcmb.20.2.3177. PMID: 9922220.
- **Satoh et al 1988:** Satoh K, Nukiwa T, Brantly M, Garver RI Jr, Hofker M, Courtney M, Crystal RG. Emphysema associated with complete absence of alpha 1- antitrypsin in serum and the homozygous inheritance [corrected] of a stop codon in an alpha 1-antitrypsin-coding exon. *Am J Hum Genet.* 1988. Jan;42(1):77-83.
- **Schaefer et al 2018:** Schaefer B, Mandorfer M, Viveiros A, Finkenstedt A, Ferenci P, Schneeberger S, Tilg H, Zoller H. Heterozygosity for the alpha-1-antitrypsin Z allele in cirrhosis is associated with more advanced disease. *Liver Transpl.* 2018;24(6):744-751. doi:10.1002/lt.25057
- **Schneider et al 2020:** Schneider CV, Hamesch K, Gross A, Mandorfer M, Moeller LS, Pereira V, Pons M, Kuca P, Reichert MC, Benini F, Burbaum B, Voss J, Gutberlet M, Woditsch V, Lindhauer C, Fromme M, Küppers J, Bewersdorf L, Schaefer B, Eslam M, Bals R, Janciauskienė S, Carvão J, Neureiter D, Zhou B, Wöran K, Bantel H, Geier A, Dirrichs T, Stickel F, Teumer A, Verbeek J, Nevens F, Govaere O, Krawczyk M, Roskams T, Haybaeck J, Lurje G, Chorostowska-Wynimko J, Genesca J, Reiberger T, Lammert F, Krag A, George J, Anstee QM, Trauner M, Datz C, Gaisa NT, Denk H, Trautwein C, Aigner E, Strnad P; European Alpha-1 Liver Study Group. Liver Phenotypes of European Adults Heterozygous or Homozygous for Pi*Z Variant of AAT (Pi*MZ vs Pi*ZZ genotype) and Noncarriers. *Gastroenterology.* 2020;159(2):534-548.e11. doi:10.1053/j.gastro.2020.04.058
- **Seersholtz et al 2000:** Seersholtz N, Wilcke JT, Kok-Jensen A, Dirksen A. Risk of hospital admission for obstructive pulmonary disease in alpha(1)-antitrypsin heterozygotes of phenotype PiMZ. *Am J Respir Crit Care Med.* 2000 Jan;161(1):81-4. doi: 10.1164/ajrccm.161.1.9812131. PMID: 10619801.

- **Seri et al 1992:** Seri M, Magi B, Cellesi C, Olia PM, Renieri A, De Marchi M. Molecular characterization of the P and I variants of alpha 1-antitrypsin. *Int J Clin Lab Res.* 1992;22(2):119-21. PubMed PMID: 1504305.
- **Sinden et al 2014:** Sinden NJ, Koura F, Stockley RA. The significance of the F variant of alpha-1-antitrypsin and unique case report of a PiFF homozygote. *BMC Pulm Med.* 2014 Aug 7;14:132. doi: 10.1186/1471-2466-14-132. PubMed PMID: 25098359.
- **Sørheim et al 2010:** Sørheim IC, Bakke P, Gulsvik A, Pillai SG, Johannessen A, Gaarder PI, Campbell EJ, Agustí A, Calverley PM, Donner CF, Make BJ, Rennard SI, Vestbo J, Wouters EF, Paré PD, Levy RD, Coxson HO, Lomas DA, Hersh CP, Silverman EK. α_1 -Antitrypsin protease inhibitor MZ heterozygosity is associated with airflow obstruction in two large cohorts. *Chest.* 2010 Nov;138(5):1125-32. doi: 10.1378/chest.10-0746. Epub 2010 Jul 1. PMID: 20595457; PMCID: PMC2972629.
- **Sproule et al 1983:** Sproule BJ, Cox DW, Hsu K, Salkie ML, Herbert FA. Pulmonary function associated with the Mmalton deficient variant of alpha 1-antitrypsin. *Am Rev Respir Dis.* 1983 Feb;127(2):237-40.
- **Stoller JK et al 2006 (updated 2017):** Stoller JK, Lacbawan FL, Aboussouan LS. Alpha-1 Antitrypsin Deficiency. 2006 Oct 27 [Updated 2017 Jan 19]. In: Adam MP, Ardinger HH, Pagon RA, et al., editors. *GeneReviews® [Internet].* Seattle (WA): University of Washington, Seattle; 1993-2020. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK1519/>
- **Stoller JK et al 2006 (updated 2020):** Stoller JK, Hupertz V, Aboussouan LS. Alpha-1 Antitrypsin Deficiency. 2006 Oct 27 [updated 2020 May 21]. In: Adam MP, Ardinger HH, Pagon RA, Wallace SE, Bean LJH, Stephens K, Amemiya A, editors. *GeneReviews® [Internet].* Seattle (WA): University of Washington, Seattle; 1993–2020. PMID: 20301692. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK1519/>
- **Strnad et al 2019:** Strnad P, Buch S, Hamesch K, Fischer J, Rosendahl J, Schmelz R, Brueckner S, Brosch M, Heimes CV, Woditsch V, Scholten D, Nischalke HD, Janciauskiene S, Mandorfer M, Trauner M, Way MJ, McQuillin A, Reichert MC, Krawczyk M, Casper M, Lammert F, Braun F, von Schönfels W, Hinz S, Burmeister G, Hellerbrand C, Teufel A, Feldman A, Schattenberg JM, Bantel H, Pathil A, Demir M, Kluwe J, Boettler T, Ridinger M, Wodarz N, Soyka M, Rietschel M, Kiefer F, Weber T, Marhenke S, Vogel A, Hinrichsen H, Canbay A, Schlattjan M, Sosnowsky K, Sarrazin C, von Felden J, Geier A, Deltenre P, Sipos B, Schafmayer C, Nothnagel M, Aigner E, Datz C, Stickel F, Morgan MY, Hampe J, Berg T, Trautwein C.. Heterozygous carriage of the alpha1-antitrypsin Pi*Z variant increases the risk to develop liver cirrhosis. *Gut.* 2019;68(6):1099-1107. doi:10.1136/gutjnl-2018-316228
- **Suh-Lailam et al 2014:** Suh-Lailam BB, Procter M, Krautscheid P, Haas J, Kumar S, Mao R, Grenache DG. Challenging identification of a novel PiSF and the rare PiMmaltonZ α_1 -antitrypsin deficiency variants in two patients. *Am J Clin Pathol.* 2014 May;141(5):742-6. doi: 10.1309/AJCPR7EIQS8PIMLV. PubMed PMID: 24713750.
- **Takabe et al 1992:** Takabe K, Seyama K, Shinada H, Nouchi T, Miyahara Y, Nukiwa T, Miyake K, Tsukimoto K, Ichioka M, Marumo F. A new variant of alpha-1-antitrypsin deficiency (Siiyama) associated with pulmonary emphysema. *Intern Med.* 1992. May;31 (5):702-7.
- **Takahashi et al 1988:** Takahashi H, Nukiwa T, Satoh K, Ogushi F, Brantly M, Fells G, Stier L, Courtney M, Crystal RG. Characterization of the gene and protein of the alpha 1-antitrypsin "deficiency" allele Mprocida. *J Biol Chem.* 1988. Oct 25; 263(30):15528-34.
- **Tete-Benissan & Gbeassor 2011:** Tete-Benissan A, Gbeassor M. Alpha-1-antitrypsin phenotypes in togolese population: description of high frequency of rare allele Pi(F) in isolated ethnic group. *Pathol Biol (Paris).* 2011 Oct;59(5):269-74. doi:

- 10.1016/j.patbio.2009.10.003. Epub 2009 Nov 25.
- **Veith et al 2019:** Veith M, Klemmer A, Anton I, El Hamss R, Rapun N, Janciauskiene S, Kotke V, Herr C, Bals R, Vogelmeier CF, Greulich T. Diagnosing Alpha-1-Antitrypsin Deficiency Using A PCR/Luminescence-Based Technology. *Int J Chron Obstruct Pulmon Dis.* 2019 Nov 18;14:2535-2542.
 - **Yuasa et al 1993:** Yuasa I, Sugimoto Y, Ichinose M, Matsumoto Y, Fukumaki Y, Sasaki T, Okada K. PI**S*(iiyama), a deficiency gene of alpha 1-antitrypsin: evidence for the occurrence in western Japan. *Jpn J Hum Genet.* 1993 Jun;38(2):185-91. PubMed PMID: 8358043.
 - **Zhumagaliyeva et al 2017:** Zhumagaliyeva A, Ottaviani S, Greulich T, Gorrini M, Vogelmeier C, Karazhanova L, Nurgazina G, DeSilvestri A, Kotke V, Barzon V, Zorzetto M, Corsico A, Ferrarotti I. Case-finding for alpha1-antitrypsin deficiency in Kazakh patients with COPD. *Multidiscip Respir Med.* 2017 12:23. doi: 10.1186/s40248-017-0104-5. PubMed PMID: 29090095
 - **Zorzetto et al 2008:** Zorzetto M, Russi E, Senn O, Imboden M, Ferrarotti I, Tinelli C, Campo I, Ottaviani S, Scabini R, von Eckardstein A, Berger W, Brändli O, Rochat T, Luisetti M, Probst-Hensch N; SAPALDIA Team. SERPINA1 gene variants in individuals from the general population with reduced alpha1-antitrypsin concentrations. *Clin Chem.* 2008 Aug;54(8):1331-8. doi: 10.1373/clinchem.2007.102798. PubMed PMID: 18515255

Trademarks:

The following product names and any registered and unregistered trademarks mentioned in this document are used for identification purposes only and remain the exclusive property of their respective owners.

ORAcollect·Dx is a trademark of DNA Genotek Inc.

Luminex's xMAP is a trademark of Luminex Inc.

AlphalD™ is a trademark of Grifols Worldwide Operations Limited.

A1AT Genotyping Test is a trademark of Progenika Biopharma S.A.

Customer Support:

While AlphalD™ At Home Genetic Health Risk Service is designed for ease and convenience, you may have questions along the way about collecting your saliva sample, registering your AlphalD™ At Home Saliva Collection Kit, or accessing your AlphalD™ At Home Genetic Health Risk Service Report. Our [FAQs & Resources](#) page is a great place to start. You can also reach our Customer Support at [1-844-891-3030](tel:1-844-891-3030). If you have any questions about your AlphalD™ At Home Genetic Health Risk Service Report and what it means to you, please speak with a healthcare professional.

Revision History:

Version PI_12923D0000_02:

- Update of the 12th warning in section “Warnings and limitations of the Service”.
- Inclusion of an assay limitation in section “Tested interfering variants” related to the tested interfering variant PI*P Donauwoerth (c.1093G>A, rs143370956).
- Inclusion of an assay limitation in section “Tested interfering variants” related to the tested interfering variant PI*Q0 bolton (c.1158delC, rs764325655). Inclusion of the reference 13 in the General References section.
- Update of the summary of the AlphaID™ At Home Saliva Collection Kit user study for ORAcollect®-Dx (OCD-100.014).

Version PI_12923D0000_03:

- Inclusion of an assay limitation in section “Tested interfering variants” related to the tested interfering variant rs1019260714.

Version PI_12923D0000_04:

- Inclusion of an assay limitation in section “Tested interfering variants” related to the tested interfering variant rs199422209.

Version PI_12923D0000_05:

- Inclusion of an assay limitation in section “Tested interfering variants” related to the tested interfering variant rs372571769.

Version PI_12923D0000_06:

- Inclusion of an assay limitation in section “Tested interfering variants” related to the tested interfering variant rs201774333.